# LIQUIZYME CALCIUM

(OCPC Method)

| Code   | Product Name                         | Pack Size |
|--------|--------------------------------------|-----------|
| LS035B | Liquizyme Calcium (OCPC)             | 500 ml    |
| LS035D | Liquizyme Calcium(OCPC)              | 100 ml    |
| LS035E | 035E Liquizyme Calcium(OCPC) 1000 ml |           |
| LS035F | Liquizyme Calcium(OCPC)              | 5 LTR     |

#### Intended Use

Diagnostic reagent for quantitative *in vitro* determination of Calcium in human serum.

#### Clinical Significance

Calcium has numerous function within the body, not only as a structural factor in bones and teeth, but also in normal neuromuscular function and the clotting of blood.

Hypercalcaemia may develop in patients with Paget's disease of bone and hyperparathyroidism. The cause of hypercalcaemia in malignancy is an increased bone resorption either caused by metastasis or by humoral factors produced by the tumor cell.

In Rickets, Coeliac diseases, idiopathic steatorrhea, osteomalacia, tropical sprue and following surgical resection of the small intestine, serum calcium is often moderately reduced, usually in association with low plasma protein concentration.

## Principle

OCPC reacts with calcium in alkaline solution to form a purple coloured complex. The intensity of the purple color formed is proportional to the concentration of calcium and is measured photometrically between 540 nm and 600 nm with maximum absorbance at 578 nm.

#### **Reagent Composition**

 $\begin{tabular}{lll} \textbf{Reagent 1: Calcium OCPC} & \textbf{Reagent} \\ \textbf{OCPC} & : & <2 \, mmol/L \\ \textbf{8-Hydroxy quinoline} & : & <1 \, mmol/L \\ \end{tabular}$ 

Reagent 2: DEA Buffer: QS

Reagent 3: Calcium Standard : 10 mg/dl

Ready to use

#### Reagent Preparation

Bring all the reagents to room temperature.

The working reagent is prepared by mixing equal parts of reagent 1 and reagent 2 (1:1 ratio).

The working reagent is ready for the use.

#### Materials Required But Not Provided

- Clean & Dry container.
- Laboratory Glass Pipettes or Micropipettes & Tips
- Colorimeter or Bio-Chemistry Analyzer.

## Stability And Storage

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at  $2-8^{\circ}$ C.



#### Specimen Collection And Handling

Use unheamolytic serum. It is recommended to follow NCCLS procedures (or similar standardized procedure).

#### Stability In Serum:

7 days : at 2-8°C

#### Calibration

Calibration with the Calcium standard provided in the kit is recommended.

#### **Quality Control**

It's recommended to run normal and abnormal control sera to validate reagent performance.

#### **Unit Conversion**

 $mg/dl \times 0.25 = mmol/L$ 

#### **Expected Values**

Serum : 9.0-10.6 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

#### Performance Data

Data contained within this section is representative of performance on Beacon system. Data obtained in your laboratory may differ from these values.

#### Precision

| Intra-assay precision | Mean    | SD      | CV   |
|-----------------------|---------|---------|------|
| Within run (n=20)     | (mg/dl) | (mg/dl) | (%)  |
| Sample 1              | 8.59    | 0.05    | 0.53 |
| Sample 2              | 11.87   | 0.15    | 1.28 |
| Inter-assay precision | Mean    | SD      | CV   |
| Run to run (n=20)     | (mg/dl) | (mg/dl) | (%)  |
| Sample 1              | 9.51    | 0.03    | 0.35 |

#### Comparison

A comparison between Beacon Calcium (y) and a commercially available test (x) using 20 samples gave following results:

y = 1.008 x - 0.084

r = 0.994

BEACON DIAGNOSTICS PVT. LTD. 424, NEW GIDC, KABILPORE, NAVSARI - 396 424. INDIA

#### Interferences

 $Following \, substances \, do \, not \, interfere \, : \,$ 

haemoglobin up to 10 g/l, bilirubin up to 40mg/dl, trigly cerides up to 500 mg/dl.

# Warning And Precautions

For  $\ensuremath{\textit{in vitro}}$  diagnostic use. To be handled by entitled and professionally educated person. Reagent of kit not classified as dangerous.

#### Waste Management

 ${\it Please \, refer \, to \, local \, legal \, requirements.}$ 

# Assay Procedure

Wavelength : 578 nm : 1 cm Cuvette

| Addition Sequence | Reagent Blank | Standard | Sample  |
|-------------------|---------------|----------|---------|
| Reagent 1         | 1000 μΙ       | 1000 μΙ  | 1000 μΙ |
| Standard          | -             | 20 μΙ    | -       |
| Sample            | -             | -        | 20 μΙ   |
| Distilled Water   | 20 μΙ         | -        | -       |

Mix well incubate at RT for 5 mins. Measure the absorbance of standard(Abs.S) And test (Abs. T) Against reagent blank at 578 nm.

# Calculation

Calcium (mg/dl) = — x 10 Abs. S

# Assay Parameters For Photometers

| Mode                   | End point        |  |
|------------------------|------------------|--|
| Wavelength 1 (nm)      | 578              |  |
| Sample Volume (μl)     | 20               |  |
| Reagent Volume (μΙ)    | 1000             |  |
| Incubation time (min.) | 5                |  |
| Incubation temp. (°C)  | Room Temperature |  |
| Normal Low (mg/dl)     | 9.0              |  |
| Normal High (mg/dl)    | 10.6             |  |
| Linearity Low (mg/dl)  | 0.6              |  |
| Linearity High (mg/dl) | 15               |  |
| Standard Concentration | 10 mg/dl         |  |
| Blank with             | Reagent          |  |
| Unit                   | mg/dl            |  |

# References

1. Bishop, M.C DUBEB-VON LAUFEN, J,L.,Burtis ,carl Aa and ashwood ,Tiet 110,61

# Symbols Used On Labels

REF

Catalogue Number

Manufacturer

Lot Number

li

See Instruction LOT for Use

CONT

Content

Storage Temperature



**Expiry Date** 



In Vitro Diagnostics





BEA/24/CAO/LS/IFU-03 20/10/2022